

FitMao Body Composition Analyzer Eight-electrode bioelectrical impedance measurement technology

BIA

**BIOELECTRICAL
IMPEDANCE**

The electrodes are applied to different parts of the human body to be measured, and the eight-electrode multi-frequency impedance acquisition method is used to measure the resistance-capacitance model of the biological ternary model. The error between the measured value and the theoretical value is extremely small, and the accuracy of the measurement system is optimal.



Why choose us?

Professional health measurement and fitness data management capabilities

	FitMao	other brands
Measurement technology	High-tech certification	common technology
Electrode	Eight electrodes	Four-electrode/ eight-electrode
Frequency	Multi-frequency five-segment	single frequency/ multi frequency
Measurement and Analysis	Direct measurement + big data analysis	Calculation of empirical variables
System Support	FitMao fitness management system	no system support
Scenes to be used	Fitness studio/club/chain store/ beauty body, etc. to improve the grade of the venue	general scene



9 characteristics

Intelligent interconnection,
precise and stable

BIA

Bioelectrical
impedance



Eight
electrodes



Multi-Frequency
Segment Test



10 inch IPS
touch screen



Multi-terminal
physical test report



Directly get
the report



personal
health record



Data cloud
service



High-tech
certification



Face recognition user information, making physical examination more convenient

In addition to the traditional login method, FM-280 Face recognition height measuring instrument fully supports the face recognition test. There is no need to manually input user information, and information matching can be completed by swiping the face. Fast and efficient.




Fitness membership card (IC card) identification

FM-280 Face recognition height measuring instrument has a built-in IC card reader module. After binding the IC card, the user can perform identity matching through the traditional gym membership card, and the physical test data will be permanently stored in the cloud.




Human Body Segmented Electrical Impedance Analysis Technology

Direct piecewise multi-frequency bioelectrical impedance analysis, treating the human body as consisting of five cylinders: left arm, right arm, left leg, right leg, and torso. FitMao support independent measurements of individual cylinders, giving accurate measurements of the whole body.




Multi-frequency (high + low frequency) direct current measurement

Different high and low frequency conditions are used to measure intracellular and extracellular fluids for accurate total body water analysis. The use of multiple frequencies improves the measurement accuracy of the body composition analyzer.



Want to track real changes in your body?
Is the body getting healthy? How has it changed?

The body composition curve will tell you!



Through the body composition analyzer, not only can we see the changes in weight gain, but also know the changes in body fat percentage over time. At the same time, it can also compare whether important physical information such as skeletal muscle weight has changed.

IPS HD touch capacitive screen

Equipped with a 10-inch high-definition screen, there is no need to install additional mobile APPs. It integrates touch operation and management, which can realize body composition analysis, view historical physical examination records, print paper reports, etc.



More than 40 human health data

Professional-grade body composition analysis report

Body Composition Analysis Report

Name 周黛墨 Height 161.3cm Age 25 Gender female Time/Date 2020/06/15 14:46

Body Composition Analysis

	Values	Body Water	Soft Lean Mass	Fat Free Mass	Weight
ECF Extracellular fluid	L 10.6 (9.5~11.8)	28.5 (26.3~32.1)	34.8 (33.8~43.1)	36.2 (35.8~43.7)	58.2 (43.9~59.5)
CF Intracellular fluid	L 15.9 (13.5~18.4)				
Protein	kg 7.4 (7.0~8.6)				
Minerals	kg 2.72 (2.44~2.98)				
Fat	kg 20.2 (10.3~16.5)				

Muscle-Fat Analysis

	Under	Normal	Over
Weight kg	30.0 35.0 40.0 43.9 45.0 50.0 55.0 59.5 65.0 70.0 75.0	58.2	
SMM Skeletal Muscle Mass %	5.0 10.0 15.0 20.0 25.0 30.0 35.0 40.0 45.0 50.0 55.0	19.6	
Body Fat Mass kg	6.0 8.0 10.0 12.0 14.0 16.0 18.0 19.0 22.0 25.0 28.0	21.8	

Obesity Analysis

	Under	Normal	Over
BMI Body Mass Index kg/m2	30.0 35.0 40.0 43.9 45.0 50.0 55.0 59.5 65.0 70.0 75.0	22.6	
PBF Percent Body Fat %	5.0 10.0 15.0 20.0 25.0 30.0 35.0 40.0 45.0 50.0 55.0	23.2	

Segmental Assessment

Segmental Lean Analysis	RA		LA	
	2.99kg	99.2%	2.97kg	98.3%
	Normal	Normal	Normal	Normal
	24.19kg	98.8%	24.19kg	98.8%
Segmental Fat Analysis	RA		LA	
	0.65kg	103.2%	0.66kg	105.2%
	Normal	Normal	Normal	Normal
	2.99kg	99.2%	2.99kg	105.2%
Segmental Lean Analysis	RL		LL	
	9.91kg	112.9%	9.68kg	110.2%
	Under	Under	Under	Under
	24.19kg	98.8%	24.19kg	98.8%
Segmental Fat Analysis	RL		LL	
	1.99kg	102.2%	1.89kg	107.1%
	Normal	Normal	Normal	Normal
	2.99kg	99.2%	2.99kg	105.2%

Body Composition Change

	20/06/02 17:31	20/06/03 15:21	20/06/04 16:25	20/06/05 17:56	20/06/06 10:36	20/06/07 11:56	20/06/09 16:17
Weight kg	62.3	61.0	60.3	59.8	59.5	58.6	58.2
SMM Skeletal Muscle Mass %	18.2	18.6	18.9	19.3	19.1	19.2	19.6
PBF Percent Body Fat %	26.3	25.1	24.6	23.9	23.1	22.5	21.8


FiTMao-280

Score 62.7 Physiological Age 26

Comprehensive Evaluation

Basal Metabolic Rate	1639	kcal
Bone mineral content	2.8	kg
Waist hip ratio	0.88	(0.79~0.90)
Visceral Fat Level	11.3	(4.0~8.0)

Body Type



Weight Control

Target Weight	51.8kg
Weight Control	- 7.4kg
Fat Control	- 9.8kg
Muscle Control	+ 2.6kg

Results QR Code



Impedance

	RA	LA	TR	RL	LL
5KHZ	325.0	327.5	27.9	241.1	255.3
50KHZ	325.0	327.5	27.9	241.1	255.3

Built-in wireless printing program

WIFI printer, one-click can print out the measurement result




Compatible with mainstream WIFI printers on the market, EPSON, HP, etc. Please purchase according to the specified brand or model, please contact customer service for details.



No need to specify special paper

The template and data are output at one time, and the overall printing can be realized, and the data is not misplaced. There is no need to specify special printing paper, which greatly saves operating costs.



Some body fat analyzer merchants will sell special printing paper



The price of the special paper is high



FitMao only need ordinary blank A4 printing paper

FitMao Product Parameters

Bioresistance Measurement Item	Fat weight, total water, protein, inorganic salts, body fat percentage, body water percentage, skeletal muscle percentage, visceral fat level, body mass index BMI, basal metabolism, bone mineral content, extracellular fluid, intracellular fluid, impedance, body weight, Lean body mass, composite score, body age, segmental fat mass (limbs + trunk), segmental muscle mass (limbs + trunk), segmental fat ratio (limbs + trunk)
Measurement methods	Direct segmental multi–frequency bioelectrical impedance analysis SDM–BIA measurement method"
Electrode method	Eight electrodes
Measurement and Analysis	Direct measurement + big data analysis
Display screen	10.1 inch HD touch capacitive screen
Report Type	Android console report,PC report,Mobile application H5 report,Paper report in A4 format
Voice broadcast function	Full live voice prompts
Data storage	Cloud data storage (need to be connected to the Internet), support PC terminal; FitMao smart front desk, fitness assistant, etc.